



SOP Title: Agilent DNA 7500 Quality Control	Version 1.1, Page 1
Authorized Signature:	Issuing Date: 02.09.11
	Last Revision Date: 11.30.11 - SV
Staff Able to Perform Procedure: Research Technician and higher	
Principle of the Method: This protocol describes the initial quality and quantification control steps when total DNA is received by the UB Next-Generation Sequencing and expression Analysis Core.	
Sample Type: total DNA sample	
Equipment Requirements: <ul style="list-style-type: none"> o Agilent 2100 Bioanalyzer (Agilent) o Bioanalyzer Chip Vortex (IKA) o Vortex (VWR) o Microcentrifuge (VWR) o Timer (1 hour) 	
Reagents & Material Requirements: <ul style="list-style-type: none"> o Quant-IT Picogreen Assay (Invitrogen) o PMMA Cuvettes (Fischer) o RNase- free centrifuge tubes 1.5ml, 200 ul (USA Scientific) o Pipette Man 2ul, 20ul, 200ul, 1000ul (Rainin) o Filter Pipette Tips (Rainin) 	
Sample <ul style="list-style-type: none"> o Wearing gloves, place total DNA samples at -80°C in labeled 5x5 Cryo freezer box. 	
Preparing the Gel-Dye Mix <ul style="list-style-type: none"> o Allow DNA dye concentrate (blue) and DNA gel matrix (red) to equilibrate to room temperature for 30 minutes. o Vortex DNA dye concentrate (blue) and add 25 µl of the dye to a DNA gel matrix vial (red). o Vortex solution well and spin down. Transfer to spin filter. o Centrifuge at 1500 g for 10 minutes. Protect solution from light. Store at 4°C. 	
Loading the Gel-Dye Mix <ul style="list-style-type: none"> o Adjust the syringe clip so it is at the top position. o Allow the gel-dye mix equilibrate to room temperature for 30 minutes before use. o Put a new DNA chip on the chip priming station. o Pipette 9.0 µl of gel-dye mix in the well marked  o Make sure that the plunger is positioned at 1 ml and then close the chip priming station. o Press plunger until it is held by the clip. o Wait exactly 30 seconds then release clip. o Wait for 5 seconds. Slowly pull back plunger to 1 ml position. o Open the chip priming station and pipette 9.0 µl of gel-dye mix in the wells marked  	
Loading the Markers <ul style="list-style-type: none"> o Pipette 5 µl of marker (green) in all 12 sample and ladder wells. Do not leave any wells empty. 	
Loading the Ladder and the Samples <ul style="list-style-type: none"> o Pipette 1 µl of DNA ladder (yellow) in the well marked ladder. o In each of the 12 sample wells pipette 1 µl of sample (used wells) or 1 µl of de-ionized water (unused wells). o Put the chip horizontally in the adapter and vortex for 1 minute at the indicated setting (2400 rpm). o Run the chip in the Agilent 2100 bioanalyzer within 5 minutes. 	