

## <u>Standard Operating Procedure</u> for FREEZING INDUCED PLURIPOTENT STEM CELLS

## 1.0 <u>PURPOSE</u>

This procedure describes how to freeze induced pluripotent stem cells (iPSC).

2.0 <u>SCOPE</u> (Should include which Cores this SOP applies to) This procedure applies to all personnel who wish to freeze iPSC.

## 3.0 <u>PROCEDURE</u>

- 1. Prepare fresh 2× freezing medium (60% FBS, 20% DMSO, 20% hES medium) and place it on ice.
- 2. Aspirate the culture medium, wash each well of the 6-well plate with 1 mL DMEM/F12.
- 3. Add 0.5 mL dispase solution to each well of 6-well plate.
- 4. Incubate the dish for 3-5 minutes in 37°C incubator.
- 5. Aspirate the dispase solution and wash each well with 2 ml DMEM/F12.
- 6. Add 2 mL of hES medium and dislodge the cells using a glass pipet. Transfer the cells to a 15 mL tube. Wash the plate with another 1 ml hES medium and transfer it to the same tube.
- 7. Centrifuge the cells at 200 g for 3 minutes at room temperature.
- 8. Discard the supernatant, resuspend the cells in hES medium (usually, 1 ml for one well). Add an equal volume of  $2\times$  freezing medium to the resuspend cell and mix them gently.
- 9. Aliquot 1 mL of the cell suspension into each cryovial. Put the cryovials into a cryo freezing container and transfer them to -80°C overnight.
- 10. In the next day, transfer the cells to a liquid nitrogen tank for long term storage.

Jian Feng 8/23/12