



Standard Operating Procedure for *Subculture IPS Cells*

1.0 PURPOSE

This procedure describes how to subculture IPS cells.

2.0 SCOPE (Should include which Cores this SOP applies to)

This procedure applies to all personnel who wish to subculture IPS cells.

3.0 PROCEDURE

1. Prepare MEF culture 6-well dishes one day before subculture iPSC.
2. Before split the iPSC, wash each well of MEF culture dishes with 2 ml DMEM/F12 medium. Add 1.5 ml hES medium containing 10 μ M Y-27632/Rock inhibitor into each well.
3. Split iPSC.
 - a. Aspirate hES medium from iPSC culture dishes. Wash each well of 6-well dishes with 1 ml DMEM/F12 medium.
 - b. Add 0.5ml dispase into each well and incubate in 37 °C incubator for 3-5 mins.
 - c. Wash each well with 2 ml DMEM/F12 medium.
 - d. Add 2 ml hES medium containing 10 μ M Y-27632/Rock inhibitor into each well.
 - e. Detach iPSC colonies from the dishes with a glass pipette and put the cells into a 15-ml tube. Rinse the well with 1 ml hES medium containing 10 μ M Y-27632/Rock inhibitor and transfer it into the same tube. Gently cut colonies into small pieces with a pipette.
 - e. Add a suitable amount of iPSC cells into each well and make sure that the total volume of each well is not less than 2 ml. Note: usually one well of iPSC colonies is split to six wells.

Jian Feng 8/23/12